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

Operations Guide of Sepaths

1 Preparation work	<p>1.1 Please Check power and PC COM lines. Then push the power button and open the software of Sepaths.</p> <p>1.2 PC-COM Connection:[System sets]-[PC COM], Please check the instrument and select the COM port number. (Tips: Please click the lamp picture .The light can be controlled and PC COM connection is successful. If not, please close the software and turn off the instrument. Please check the lines and repeat above operations) .</p> <p>1.3 Check Solvents:[System sets]-[Solvent set],Please select the solvent names for eight solvent bottles and make sure the volume of solvent is enough to the experiment.</p> <p>1.4 Set up pressure value: (Please adjust the pressure from lower pressure to the aim- pressure)</p> <p>1) Input pressure: Please adjust the safe valve and make sure the input pressure is below 60PSI/0.4Mpa. Make sure the gas pipeline without gas leak .</p> <p>2) Solvent pressure: About 15 PSI/0.1Mpa;</p> <p>3) Extraction pressure (Below pressure value is from drinking water and for reference only, can be adjusted.)</p> <p>A) Sample volume<=20ml /Cartridge or disk /About 10PSI</p> <p>B) Sample volume >20ml /Cartridge only/About 4-15PSI</p> <p>C) Sample volume >20ml /Disk only/About 10-15PSI</p> <p>1.5 Prepare sample: If Volume<=20ml, please install small volume adapter. If Volume>25ml, please use sample bottle only.</p> <p>1.6 Prepare vial of collection: Please wash and dry the vials. There are two collection positions. Position one is fixed and position two can be choosed .</p> <p>1.7 Prepare cartridges and disks: Please wash cartridge adapters and disk adapters with cleaning water and ethanol before and after every experiment.</p> <p>1.8 Check waste liquid barrel: There are three waste liquid tubes. First is waste water tube. Another two are waste organic solvents tubes.</p>
2 New and Load Method	<p>2.1 New Method: Please choose [Extraction Type] first. The sample volume must be quantified, because Sepaths only treat the whole sample.</p> <p>2.2 Load Method: If we do not have the method, please new a method. If we have a method, please load the method directly. Each channel can choose the same or different method. And from the channel one to channel six, the same method will be run at the same time. For example, channel one and channel three are both method two, and channel two and channel four are both method one. After click start button, the method two will be run first, the method one will be run second.</p> <p>2.3 Run Method: After loading method and everything is ok, Please click start button and the method will be run. If we click the stop button, we can choose stop method or pause method. If we have chosen pause method, click the start button and the method will go on. When every step of SPE method is completed, Sepaths will tell the operator through buzzer.</p> <p>2.4 Check sampling step: We can choose to pause the method after the step of sampling automatically and to check the sample bottles. If sample bottles are empty, please click button of Sure and the next steps will be going on.</p> <p>2.5 Chinese and English software surface: It is important that English methods should only be used in English software and Chinese method should be used in English software.</p>
3 Finish work	<p>3.1 Check the volume of collected liquid and make sure it is right.</p> <p>3.2 Please close the solvent pressure and extractor pressure .If it is necessary, we need close the valve of nitrogen.</p> <p>3.3 Please uninstall the adpters of cartridges or disks and wash with water and ethanol in order to avoid the cross contamination.</p> <p>3.4 Please turn off the software and instrument.</p> <p>3.5 Finally, clean the lab table and leave with smile.</p>

Tips: If you have any questions about Sepaths, phone number 138-1069-2460 (Mr Xie) or email xgxie@labtechgroup.cn